bs-4700R

[Primary Antibody]

Bioss ANTIBODIES

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

EBV LMP2A Rabbit pAb

- DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Target: EBV LMP2A

Immunogen: KLH conjugated synthetic peptide derived from HHV4t2 LMP2A:

401-497/497.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)

Reactivity: (predicted: HHV4t2)

Predicted MW.: 53 kDa

- SELECTED CITATIONS -

• [IF=3.329] Hui D et al. CD44+CD24-/low sphere-forming cells of EBV-associated gastric carcinomas show immunosuppressive effects and induce Tregs partially through production of PGE2. Exp Cell Res. 2020 May 15;390(2):111968. WB;human. 32197932