

bs-4700R

[Primary Antibody]

EBV LMP2A Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:50-200) ELISA (1:5000-10000) Reactivity: (predicted: HHV4t2) Predicted MW.: 53 kDa
Clonality: Polyclonal		
Target: EBV LMP2A		
Immunogen: KLH conjugated synthetic peptide derived from HHV4t2 LMP2A: 401-497/497.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=3.329]** Hui D et al. CD44+CD24-/low sphere-forming cells of EBV-associated gastric carcinomas show immunosuppressive effects and induce Tregs partially through production of PGE2. Exp Cell Res. 2020 May 15;390(2):111968. WB ;human. 32197932