

bs-11053R

[Primary Antibody]

KAL1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Rat, Pig, Sheep, Chicken, Dog, Horse) Predicted MW.: 74 kDa Subcellular Location: Secreted ,Extracellular matrix ,Cell membrane
Clonality: Polyclonal		
GeneID: 3730	SWISS: P23352	
Target: KAL1		
Immunogen: KLH conjugated synthetic peptide derived from human AF-6/l Afadin: 131-230/680.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: May be an adhesion-like molecule with anti-protease activity.		