

bs-19588R

[Primary Antibody]

PTPLAD2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: Human (predicted: Mouse, Rat, Pig) Predicted MW.: 28 kDa Subcellular Location: Cytoplasm
Clonality: Polyclonal		
GeneID: 401494	SWISS: Q5VWC8	
Target: PTPLAD2		
Immunogen: KLH conjugated synthetic peptide derived from human PTPLAD2: 1-100/232. < Cytoplasmic >		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: PTPLAD2 is a protein that belongs to the PTPLA family. The function of PTPLAD2 remains unknown.		

— SELECTED CITATIONS —

- **[IF=7.7]** Peiran Sha. et al. Ablation of LKB1 gene changes the lipid profiles of goat intramuscular fat and enhances polyunsaturated fatty acids deposition. INT J BIOL MACROMOL. 2025 Jan;285:138268 **WB ;Goat**. 39631610