

bs-18607R

[Primary Antibody]

MPEG1/Macrophage Specific Gene Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat) Predicted MW.: 76 kDa Subcellular Location: Cell membrane
Clonality: Polyclonal		
GeneID: 219972	SWISS: Q2M385	
Target: MPEG1/Macrophage Specific Gene		
Immunogen: KLH conjugated synthetic peptide derived from human MPEG1/Macrophage Specific Gene: 101-200/716. < Extracellular >		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		