

---

## ELOVL6 Rabbit pAb

Catalog Number: bs-7055R

Target Protein: ELOVL6

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Pig, Cow, Chicken)

Predicted MW: 31 kDa

Entrez Gene: 79071

Swiss Prot: Q9H5J4

Source: KLH conjugated synthetic peptide derived from human ELOVL6: 21-120/265.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Fatty acid elongase specific to C12-C16 saturated and monoinsaturated fatty acids.