

bs-7055R

[Primary Antibody]

ELOVL6 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 79071</p> <p>Target: ELOVL6</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human ELOVL6: 21-120/265.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: Fatty acid elongase specific to C12-C16 saturated and monoinsaturated fatty acids.</p>	<p>Isotype: IgG</p> <p>SWISS: Q9H5J4</p>	<p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Cow, Chicken)</p> <p>Predicted MW.: 31 kDa</p> <p>Subcellular Location: Cell membrane ,Cytoplasm</p>
---	--	---