
FVT1 Rabbit pAb

Catalog Number: bs-13233R

Target Protein: FVT1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog)

Predicted MW: 33 kDa

Subcellular: Cell membrane ,Cytoplasm

Locations:

Entrez Gene: 2531

Swiss Prot: Q06136

Source: KLH conjugated synthetic peptide derived from human FVT1: 41-140/332.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Catalyzes the reduction of 3-ketodihydrosphingosine (KDS) to dihydrosphingosine (DHS).