

bs-18481R

[Primary Antibody]

ZDHHC22/C14orf59 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 283576</p> <p>Target: ZDHHC22/C14orf59</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human ZDHHC22/C14orf59: 181-263/263.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p>	<p>Isotype: IgG</p> <p>SWISS: Q8N966</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog, Horse)</p> <p>Predicted MW.: 29 kDa</p> <p>Subcellular Location: Cell membrane</p>
--	--	---