

bs-19133R

[Primary Antibody]

IL12B Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: IL12B

Immunogen: KLH conjugated synthetic peptide derived from human IL12B: 201-300/328.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Applications: **WB** (1:500-2000)

ELISA (1:5000-10000)

Reactivity: Rat

(predicted: Human,
Mouse)

Predicted MW.: 46 kDa