

bs-6611R

[Primary Antibody]

PLRP1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 5407</p> <p>Target: PLRP1</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human PLRP1/PNLIPRP1: 375-467/467.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: PNLIPRP1 belongs to the AB hydrolase superfamily, Lipase family and contains one PLAT domain. PNLIPRP1 has a catalytic activity of: Triacylglycerol + H₂O = diacylglycerol + a carboxylate. It is expressed in the pancreas and there are three named isoforms.</p>	<p>Isotype: IgG</p> <p>SWISS: P54315</p> <p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog, Horse)</p> <p>Predicted MW.: 50 kDa</p> <p>Subcellular Location: Secreted</p>
--	---