

bs-11788R

[Primary Antibody]

PRNP Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

| | | |
|--|--|--|
| <p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 5621</p> <p>Target: PRNP</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human PRNP: 23-120/253.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: The function of PrP is still under debate. May play a role in neuronal development and synaptic plasticity. May be required for neuronal myelin sheath maintenance. May play a role in iron uptake and iron homeostasis (By similarity). Isoform 2 may act as a growth suppressor by arresting the cell cycle at the G0/G1 phase. Soluble oligomers are toxic to cultured neuroblastoma cells and induce apoptosis (in vitro).</p> | <p>Isotype: IgG</p> <p>SWISS: P04156</p> | <p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Cow, Horse)</p> <p>Predicted MW.: 25 kDa</p> <p>Subcellular Location: Cell membrane ,Cytoplasm ,Nucleus</p> |
|--|--|--|