

bs-18154R**[Primary Antibody]****phospho-IL3RB (Tyr593) Rabbit pAb**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ELISA (1:5000-10000)
Clonality: Polyclonal		
GeneID: 1439	SWISS: P32927	Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Horse)
Target: IL3RB (Tyr593)		
Immunogen: KLH conjugated synthesised acetylpeptide derived from human IL3RB around the acetylation site of Tyr593: GP(p-Y)LG.		
Purification: affinity purified by Protein A		Predicted MW.: 97 kDa
Concentration: 1mg/ml		Subcellular Location: Cell membrane
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: GM-CSF Receptor Beta is the v subunit of the heterodimeric receptor for colony stimulating factor 2, a cytokine which controls the production, differentiation, and function of granulocytes and macrophages. The encoded protein is a member of the cytokine family of receptors. This gene is found in the pseudoautosomal region (PAR) of the X and Y chromosomes. Six transcript variants encoding five different isoforms have been found for this gene, with some of the isoforms being membrane-bound and others being soluble.		