

bs-15443R

[Primary Antibody]

HEATR8 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: ELISA (1:5000-10000)
Clonality: Polyclonal		Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Dog)
GeneID: 374977	SWISS: Q68CQ1	
Target: HEATR8		Predicted MW.: 146 kDa
Immunogen: KLH conjugated synthetic peptide derived from human HEATR8: 451-550/1324.		Subcellular Location: Cell membrane
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: MROH7 (Maestro Heat Like Repeat Family Member 7) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include binding. An important paralog of this gene is MROH7-TTC4.		