

---

## DAP3 Rabbit pAb

Catalog Number: bs-7625R

Target Protein: DAP3

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit)

Predicted MW: 42 kDa

Entrez Gene: 7818

Source: KLH conjugated synthetic peptide derived from human DAP3: 151-250/398.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Involved in mediating interferon-gamma-induced cell death.