

bs-6608R

[Primary Antibody]

DDCT Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Cow, Horse) Predicted MW.: 13 kDa Subcellular Location: Cytoplasm
Clonality: Polyclonal		
GeneID: 1652	SWISS: P30046	
Target: DDCT		
Immunogen: KLH conjugated synthetic peptide derived from human DDCT: 51-118/118.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: Tautomerization of D-dopachrome with decarboxylation to give 5,6-dihydroxyindole (DHI).		