

bs-20400R

[Primary Antibody]

ZNF749 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) Reactivity: (predicted: Human) Predicted MW.: 90 kDa Subcellular Location: Nucleus
Clonality: Polyclonal		
GeneID: 388567	SWISS: O43361	
Target: ZNF749		
Immunogen: KLH conjugated synthetic peptide derived from human ZNF749: 151-250/778.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: May be involved in transcriptional regulation.		