

---

## MEI-1 Rabbit pAb

Catalog Number: bs-6387R

Target Protein: MEI-1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Pig, Sheep, Cow, Dog, Horse)

Predicted MW: 141 kDa

Entrez Gene: 150365

Swiss Prot: Q5TIA1

Source: KLH conjugated synthetic peptide derived from human MEI-1: 266-370/1274.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: MEI1 is required for normal meiotic chromosome synapsis. It may be involved in the formation of meiotic double-strand breaks (DSBs) in spermatocytes. There are 7 known isoforms.