

**bs-6387R**

**[ Primary Antibody ]**

## MEI-1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> <b>WB</b> (1:500-2000) <b>ELISA</b> (1:5000-10000)
<b>Clonality:</b> Polyclonal		
<b>GeneID:</b> 150365	<b>SWISS:</b> Q5TIA1	<b>Reactivity:</b> (predicted: Human, Mouse, Rat, Pig, Sheep, Cow, Dog, Horse)
<b>Target:</b> MEI-1		
<b>Immunogen:</b> KLH conjugated synthetic peptide derived from human MEI-1: 266-370/1274.		
<b>Purification:</b> affinity purified by Protein A		<b>Predicted MW.:</b> 141 kDa
<b>Concentration:</b> 1mg/ml		<b>Subcellular Location:</b> Cytoplasm ,Nucleus
<b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
<b>Background:</b> MEI1 is required for normal meiotic chromosome synapsis. It may be involved in the formation of meiotic double-strand breaks (DSBs) in spermatocytes. There are 7 known isoforms.		