

bs-6387R

[Primary Antibody]

MEI-1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ELISA (1:5000-10000)
Clonality: Polyclonal		
GeneID: 150365	SWISS: Q5TIA1	Reactivity: (predicted: Human, Mouse, Rat, Pig, Sheep, Cow, Dog, Horse)
Target: MEI-1		
Immunogen: KLH conjugated synthetic peptide derived from human MEI-1: 266-370/1274.		
Purification: affinity purified by Protein A		Predicted MW.: 141 kDa
Concentration: 1mg/ml		Subcellular Location: Cytoplasm ,Nucleus
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: MEI1 is required for normal meiotic chromosome synapsis. It may be involved in the formation of meiotic double-strand breaks (DSBs) in spermatocytes. There are 7 known isoforms.		