
PTRH2 Rabbit pAb

Catalog Number: bs-7101R

Target Protein: PTRH2

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Cow, Chicken, Dog, Horse)

Predicted MW: 13 kDa

Entrez Gene: 51651

Source: KLH conjugated synthetic peptide derived from human PTRH2: 71-160/179.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: The natural substrate for this enzyme may be peptidyl-tRNAs which drop off the ribosome during protein synthesis. Promotes caspase-independent apoptosis by regulating the function of two transcriptional regulators, AES and TLE1.