

bs-7101R

[Primary Antibody]

PTRH2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 51651</p> <p>Target: PTRH2</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human PTRH2: 71-160/179.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: The natural substrate for this enzyme may be peptidyl-tRNAs which drop off the ribosome during protein synthesis. Promotes caspase-independent apoptosis by regulating the function of two transcriptional regulators, AES and TLE1.</p>	<p>Isotype: IgG</p> <p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Cow, Chicken, Dog, Horse)</p> <p>Predicted MW.: 13 kDa</p> <p>Subcellular Location: Cytoplasm</p>
---	--