
PNPT1 Rabbit pAb

Catalog Number: bs-7100R

Target Protein: PNPT1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Sheep, Cow, Dog, Horse)

Predicted MW: 81 kDa

Entrez Gene: 87178

Swiss Prot: Q8TCS8

Source: KLH conjugated synthetic peptide derived from human PNPT1: 551-650/783.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: PNPT1 is a subunit of the exosome complex, which is involved in 3-prime-to-5-prime exoribonuclease activity for RNA processing and degradation. It contains one KH domain and one S1 motif domain.