

bs-19096R

[Primary Antibody]

NECAP2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 55707</p> <p>Target: NECAP2</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human NECAP2: 181-263/263.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: This gene likely encodes a member of the adaptin-ear-binding coat-associated protein family. Studies of a similar protein in rat suggest a role in clathrin-mediated endocytosis. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2009]</p>	<p>Isotype: IgG</p> <p>SWISS: Q9NVZ3</p> <p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: Mouse (predicted: Human, Rat, Rabbit, Pig, Cow, Horse)</p> <p>Predicted MW.: 28 kDa</p> <p>Subcellular Location: Cell membrane ,Cytoplasm</p>
---	---