

**bs-0807R**

**[ Primary Antibody ]**

## **Os05g0477200 protein Rabbit pAb**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### **— DATASHEET —**

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** Os05g0477200 protein

**Immunogen:** KLH conjugated synthetic peptide derived from Oryza sativa  
Os05g0477200 protein: 221-320/472.

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%  
Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated  
freeze/thaw cycles.

**Applications:** **WB** (1:500-2000)

**IHC-P** (1:100-500)

**IHC-F** (1:100-500)

**IF** (1:100-500)

**ELISA** (1:5000-10000)

**Reactivity:** (predicted: Rice)

**Predicted MW.:** 52 kDa