
TIPIN Rabbit pAb

Catalog Number: bs-19937R

Target Protein: TIPIN

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Sheep, Cow, Dog)

Predicted MW: 35 kDa

Entrez Gene: 54962

Swiss Prot: Q9BVW5

Source: KLH conjugated synthetic peptide derived from human TIPIN: 101-200/301.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.