
ICOS Rabbit pAb

Catalog Number: bs-2583R

Target Protein: ICOS

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), ELISA (1:5000-10000)

Reactivity: Human (predicted:Mouse, Rat, Pig)

Predicted MW: 20 kDa

Detected MW: 47-57 kDa

Entrez Gene: 29851

Swiss Prot: Q9Y6W8

Source: KLH conjugated synthetic peptide derived from human ICOS: 51-150/199.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: The protein encoded by this gene belongs to the CD28 and CTLA-4 cell-surface receptor family. It forms homodimers and plays an important role in cell-cell signaling, immune responses, and regulation of cell proliferation.[provided by RefSeq, Jul 2008].