

bs-2583R

[Primary Antibody]

ICOS Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ELISA (1:5000-10000)
Clonality: Polyclonal		
GeneID: 29851	SWISS: Q9Y6W8	Reactivity: Human (predicted: Mouse, Rat, Pig)
Target: ICOS		
Immunogen: KLH conjugated synthetic peptide derived from human ICOS: 51-150/199. < Extracellular >		
Purification: affinity purified by Protein A		Predicted MW.: 20 kDa
Concentration: 1mg/ml		Subcellular Location: Secreted ,Cell membrane
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: The protein encoded by this gene belongs to the CD28 and CTLA-4 cell-surface receptor family. It forms homodimers and plays an important role in cell-cell signaling, immune responses, and regulation of cell proliferation.[provided by RefSeq, Jul 2008].		