
CT45A3 Rabbit pAb

Catalog Number: bs-14092R

Target Protein: CT45A3

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), ELISA (1:5000-10000)

Reactivity: (predicted:Human)

Predicted MW: 21 kDa

Entrez Gene: 441519

Swiss Prot: Q8NHU0

Source: KLH conjugated synthetic peptide derived from human CT45A3: 21-120/189.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.