

bs-14236R

[Primary Antibody]

DDX60L Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human) Predicted MW.: 198 kDa Subcellular Location: Cytoplasm
Clonality: Polyclonal		
GeneID: 91351	SWISS: Q5H9U9	
Target: DDX60L		
Immunogen: KLH conjugated synthetic peptide derived from human DDX60L: 501-600/1706.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=6.244]** Ye Z. et al. Associations of DDX60L With the Clinical Features and Prognosis of Hepatocellular Carcinoma.. Front Oncol. 2022 Feb;12:761021-761021 IHC ;Human. 35223465