
Goat Anti-Mouse IgG H&L, SA conjugated

Catalog Number: bs-0296G-SA

Target Protein: Goat Anti-Mouse IgG H&L

Concentration: 1.0 mg/ml

Form: Liquid

Host: Goat

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),
ELISA (1:500-1000)

Reactivity: Mouse

Purification: affinity purified by Protein G

Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol.

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.