

bs-0296G-SA

[Secondary Antibodies]

Goat Anti-Mouse IgG H&L, SA conjugated



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Goat	Isotype: IgG	Applications: WB (1:500-2000)
Clonality: Polyclonal		IHC-P (1:100-500)
Target: Goat Anti-Mouse IgG H&L		IHC-F (1:100-500)
Purification: affinity purified by Protein G		IF (1:100-500)
Concentration: 1.0 mg/ml		ICC/IF (1:100-500)
Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		ELISA (1:500-1000)
		Reactivity: Mouse