

**bs-0296G-SA**

**[ Secondary Antibodies ]**

## Goat Anti-Mouse IgG H&L, SA conjugated



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<b>Host:</b> Goat <b>Clonality:</b> Polyclonal <b>Target:</b> Goat Anti-Mouse IgG H&L <b>Purification:</b> affinity purified by Protein G <b>Concentration:</b> 1.0 mg/ml <b>Storage:</b> 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.	<b>Isotype:</b> IgG	<b>Reactivity:</b> Mouse
--	---------------------	--------------------------