bs-0296G-SA

## [ Secondary Antibodies ]

## BIOSS

## Goat Anti-Mouse IgG H&L, SA conjugated

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

Reactivity: Mouse

- DATASHEET -

Host: Goat Isotype: IgG

Clonality: Polyclonal

**Target:** Goat Anti-Mouse IgG H&L **Purification:** affinity purified by Protein G

Concentration: 1.0 mg/ml

**Storage:** 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50%

glycerol.

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## - SELECTED CITATIONS -

• [IF=3.4] Jianfei Wu. et al. A potential strategy for improving offspring behavior in maternal immune activation:

Amantadine-mediated suppression of neuroinflammation. Animal Models and Experimental Medicine. 2025 Jul;: WB

;Mouse. 40697163