

bs-8191R

[Primary Antibody]

CSTP1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal GeneID: 55313 Target: CSTP1 Immunogen: KLH conjugated synthetic peptide derived from human CSTP1: 131-314/314. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: CSMD3 belongs to the CSMD family and contains fourteen CUB domains and twenty eight Sushi (CCP/SCR) domains. It is weakly expressed in most tissues, except in brain; expressed at intermediate level in brain, including cerebellum, substantia nigra, thalamus, spinal cord, hippocampus and fetal brain. Also expressed in testis. There are five named isoforms.	Isotype: IgG SWISS: Q9BRF8 Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Chicken, Dog, Horse) Predicted MW.: 36 kDa Subcellular Location: Cytoplasm
--	---