

**bs-4802R**

**[ Primary Antibody ]**

## **Salmonella paratyphi a Rabbit pAb**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Applications:** ELISA (1:5000-10000)

**Clonality:** Polyclonal

**Reactivity:** (predicted: Mouse,  
Salmonella)

**Target:** Salmonella paratyphi a

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.