

---

## HMG20A Rabbit pAb

Catalog Number: bs-18050R

Target Protein: HMG20A

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),  
ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Pig, Sheep, Cow, Chicken, Dog, Horse, Rhesus monkey,  
Gorilla, Orangutan)

Predicted MW: 40 kDa

Subcellular Nucleus

Locations:

Entrez Gene: 10363

Swiss Prot: Q9NP66

Source: KLH conjugated synthetic peptide derived from human HMG20A: 201-300/347.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.