

bs-18050R

[Primary Antibody]

HMG20A Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Pig, Sheep, Cow, Chicken, Dog, Horse, Rhesus monkey, Gorilla, Orangutan) Predicted MW.: 40 kDa Subcellular Location: Nucleus
Clonality: Polyclonal		
GeneID: 10363	SWISS: Q9NP66	
Target: HMG20A		
Immunogen: KLH conjugated synthetic peptide derived from human HMG20A: 201-300/347.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		