

**bsm-41049M**

**[ Primary Antibody ]**

## FMDV VP1 protein (type O) Mouse mAb

**Bioss**  
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn


techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<b>Host:</b> Mouse	<b>Isotype:</b> IgG1	<b>Applications:</b> <b>WB</b> (1:500-2000) <b>ELISA</b> (1:5000-10000)  <b>Reactivity:</b> FMDV
<b>Clonality:</b> Monoclonal	<b>CloneNo.:</b> 2F11	
<b>Target:</b> FMDV VP1 protein (type O)		
<b>Purification:</b> affinity purified by Protein G		
<b>Concentration:</b> 1mg/ml		
<b>Storage:</b> Size : 50ul/100ul/200ul 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Size : 200ug (PBS only) 0.01M PBS Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

### — VALIDATION IMAGES —



Sample: Lane 1: Recombinant FMDV VP1 (type O) protein, Sumo & His  
Primary: Anti-FMDV VP1 protein (type O) (bsm-41049M) at 1/1000 dilution  
Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution  
Predicted band size: kDa  
Observed band size: 46 kDa