

bs-5799R

[Primary Antibody]

BAG3 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Chicken, Dog, Horse) Predicted MW.: 61 kDa Subcellular Location: Cytoplasm
Clonality: Polyclonal		
GeneID: 9531	SWISS: O95817	
Target: BAG3		
Immunogen: KLH conjugated synthetic peptide derived from human BAG3: 441-540/575.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release. Has anti-apoptotic activity.		

— SELECTED CITATIONS —

- **[IF=1.763]** Jie Wei. et al. The Regulation of Prototype Foamy Virus 5'Long Terminal Repeats and Internal Promoter by Endogenous Transcription Factors. Intervirology. ;:1-12 WB ;Human. 34438397