
Myoglobin Rabbit pAb

Cat. No:	W0801
Name:	Myoglobin Rabbit pAb
Storage:	Shipped at 4°C. Store at -20°C for three years. Avoid repeated freeze/thaw cycles.
Clone type:	Polyclonal
Clone No:	
Host/Source:	Rabbit
Subtype:	IgG
Reactivity:	(predicted: Human)
Format:	Liquid
Concentration:	>1mg/ml
Buffer:	0.01M PBS (pH7.4) with 0.02% Proclin300.
Applications:	ELISA
Antigen:	Purified Human Myoglobin (Natural Human Myoglobin protein) (bs-41107P)