

bsm-0990M

[Primary Antibody]

HPV16 E6 Mouse mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Mouse

Isotype: IgG

Clonality: Monoclonal

CloneNo.: 2E3

Target: HPV16 E6

Immunogen: KLH conjugated synthetic peptide derived from HPV16 E6 protein: 85-158/158.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: Size : 50ul/100ul/200ul
0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Size : 200ug (PBS only)
0.01M PBS
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ELISA (1:5000-10000)

Reactivity: (predicted: HPV16)

Predicted MW.: 11 kDa