

bs-10204R

[Primary Antibody]

HPV18 E7 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: HPV18) Predicted MW.: 12 kDa
Clonality: Polyclonal		
Target: HPV18 E7		
Immunogen: KLH conjugated synthetic peptide derived from HPV18 E7 Protein: 3-50/105.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=9.3]** Fu-Chun Huo. et al. HPV E7-driven ALKBH5 promotes cervical cancer progression by modulating m6A modification of PAK5. PHARMACOL RES. 2023 Jul;;106863 WB ;Human. 37480971