
Cystatin-C Mouse mAb

Cat. No:	V0504
Name:	Cystatin-C Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	7F9L
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted: Human)
Format:	Liquid
Concentration:	≥1mg/ml
Buffer:	0.01M PBS (pH7.4) with 0.02% Proclin300.
Applications:	ELISA
Antigen:	Recombinant human Cystatin-C (bs-41175P)