

NGAL/Lipocalin 2 Mouse mAb

Cat. No:	V0601
Name:	NGAL/Lipocalin 2 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3D10N
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Format:	Liquid
Concentration:	≥1mg/ml
Buffer:	0.01M TBS (pH7.4) with 0.02% Proclin300.
Applications:	ELISA, ICA
Antigen:	Recombinant human NGAL protein:29-198/198 (bs-43518P)
Pairing Information:	

ELISA		Detection						
		V0601	V0602	V0603	V0604	V0605	V0606	V0607
Capture	V0601	-	-	-	-	-	-	-
	V0602	-	-	-	-	-	-	-
	V0603	-	-	-	-	-	-	-
	V0604	-	+	-	-	-	-	-
	V0605	+	+	+++	-	-	-	-
	V0606	+	+	-	-	-	-	-
	V0607	-	-	-	-	-	-	-
ICA		Detection						
		V0601	V0602	V0603	V0604	V0605	V0606	V0607
Capture	V0601	-	-	-	-	-	-	-
	V0602	-	-	-	-	-	-	-
	V0603	-	-	-	-	-	-	-
	V0604	-	-	-	-	-	-	-
	V0605	-	-	-	-	-	-	-
	V0606	+++	+	+	-	-	-	-
	V0607	+	+++	+	-	-	-	-