
CIB1 Rabbit pAb

Catalog Number: bs-1682R

Target Protein: CIB1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Pig, Cow, Dog, GuineaPig, Horse)

Predicted MW: 22 kDa

Entrez Gene: 10519

Swiss Prot: Q99828

Source: KLH conjugated synthetic peptide derived from human CIB1: 101-191/191.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: CIB1(Calcium and integrin binding 1) may convert the inactive conformation of integrin alpha-IIb/beta3 to an active form through the binding to the integrin cytoplasmic domain.