

bs-2325R

[Primary Antibody]

Histamine Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Applications: ELISA (1:5000-10000)

Clonality: Polyclonal

Reactivity: (predicted: Histamine)

Target: Histamine

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Subcellular Location: Secreted

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Histamine is located in mast cells, endocrine cells of the gut, blood cells and in some cells of the peripheral and central nervous system. Histamine is a potent vasodilator when secreted by mast cells found in various tissues as a result of allergic hypersensitivity or inflammation. Histamine causes running nose, sneezing and itching, and narrows the airways in the lungs. In the central nervous system, Histamine is putative neurotransmitter. In the brain, its highest content has been found in the hypothalamus and in certain areas of the mesencephalon.
CAS: 51-45-6