

**bs-8850R**

**[ Primary Antibody ]**

## Acetyl Lysine Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<b>Host:</b> Rabbit <b>Clonality:</b> Polyclonal <b>Target:</b> Acetyl Lysine <b>Purification:</b> affinity purified by Protein A <b>Concentration:</b> 1mg/ml <b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.	<b>Isotype:</b> IgG	<b>Applications:</b> <b>IHC-P</b> (1:100-500) <b>IHC-F</b> (1:100-500) <b>IF</b> (1:100-500) <b>ICC/IF</b> (1:100-500) <b>ELISA</b> (1:5000-10000)  <b>Reactivity:</b> Species independent
--	---------------------	--

### — SELECTED CITATIONS —

- **[IF=9.417]** Zhenyin Chen. et al. Uninterrupted dynamic stiffening microenvironment enhances the paracrine function of mesenchymal stem cells for vascularization through chromatin remodeling. MATER DESIGN. 2022 Dec;224:111328 IF ;Rat. 10.1016/j.matdes.2022.111328