

bs-19119R**[Primary Antibody]****MYADML2 Rabbit pAb**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat) Predicted MW.: 33 kDa Subcellular Location: Cell membrane
Clonality: Polyclonal		
GeneID: 255275	SWISS: A6NDP7	
Target: MYADML2		
Immunogen: KLH conjugated synthetic peptide derived from human MYADML2: 101-200/307.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=7.3]** Quanyu Chen. et al. Comparative effects of hepatocyte growth factor and tacrolimus on acute liver allograft early tolerance. FRONT IMMUNOL. 2023; 14: 1162439 IF ;Rat. 37614233
- **[IF=6.107]** Bei Zhang. et al. CRISPR activation screening in a mouse model for drivers of hepatocellular carcinoma growth and metastasis. ISCIENCE. 2023 Mar;26:106099 IHC ;Human. 36843840