
LIN37 Rabbit pAb

Catalog Number: bs-18278R

Target Protein: LIN37

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog, Horse)

Predicted MW: 28 kDa

Subcellular: Cytoplasm

Locations:

Entrez Gene: 55957

Swiss Prot: Q96GY3

Source: KLH conjugated synthetic peptide derived from human LIN37: 171-246/246.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a protein expressed in the eye. [provided by RefSeq, Jul 2008]